




# tweak

## Model Comparison

	FOCUS	FOCUS+T
Available in <b>2 Colors</b>	Beige and Metallic Grey	
<b>Slim profile case, with thin tube coupling to the ear</b> comfortable, nearly invisible, no barrel effect	✓	✓
<b>Sophisticated, patented technology</b> allows you to choose your own amplification profile from 3 acoustic algorithms	✓	✓
<b>Cutting-edge automatic feedback control</b> minimizes whistling	✓	✓
 <b>Omnidirectional mode (Quiet program)</b> for all-around sound	✓	✓
 <b>Directional mode (Noise program)</b> allows for better listening in noisy environments— reduces background noise and clatter <i>(2 microphones work together to provide true directionality)</i>	✓	✓
 <b>Telecoil</b> allows you to hear better on the telephone or connect to a loop system		✓
<b>Digital volume control</b> fine tuning volume within your amplification profile for flexibility in difficult listening situations	✓	✓
<b>12 band graphic equalizer</b> —provides smooth, natural sound with amplification only where you need it	✓	✓
<b>Fast-acting sound-activated compression</b> (4 channels—improves audibility of soft sounds and softens sudden loud sounds)	✓	✓
<b>User reset capability</b>	✓	✓
<b>Low battery tone indicator</b>	✓	✓
<b>Battery size</b>	<b>312</b>	<b>312</b>
<b>Push-button access to listening modes:</b> Quiet, Noise, and Telecoil allow you to quickly and easily customize your listening preference, depending on your situation	<b>Quiet Noise</b>	<b>Quiet Noise Telecoil</b>

